

MERLIN MICROSHOT SYRINGE COMPATIBILITY GUIDE

Note: A 5 uL syringe will deliver half the nominal volume.

All syringes in this table have 42 mm needles.

Volume (uL)	Needle Ga.	Termination	Description (other)	Agilent P/N	Hamilton P/N	SGE P/N	ILS P/N
5	23	Fixed		G4513-80213 9301-0892	87987	001810	2104305
5	23	Fixed	Flexible Plunger				2164305
5	23	Removable		G4513-80234 5182-0834	87957	001815	2104309
5	23	Removable	Flexible Plunger				2164309
5	26	Fixed		G4513-80226 9301-0891	87988	001800	2104303
5	26	Removable			87958	001805	2104307
10	23	Fixed		G4513-80209 9301-0713	80387	002810	2104405
10	23	Fixed	Flexible Plunger			002835	2164405
10	23	Removable		G4513-80235 5181-8806	80357	002815	2104409
10	23	Removable	Flexible Plunger			002845	2164409
10	26	Fixed		G4513-80216 9301-0714	80388	002800	2104403
10	26	Removable			80358	002805	2104407